

## TECHNICAL SPECIFICATION

Oxygen Flow	0.5 to 5 Liters per minute
Delivery Method	Continuous
Oxygen Purity	93% ± 3%
Weight (kg)	15.2kg
Size (LxWxH)	370 x 330 x 580mm
Oxygen Outlet Pressure	8.5psi ± 0.5psi
Sound Levels (dBA)	≤40
Cord Length (m)	2 m
Electrical Requirements	220V/50Hz

# OXYGEN CONCENTRATOR

- 1-5L/min.
- Fast & easy-maintenance.
- Light-weight and robust design.
- Oxygen and performance monitor.
- Suitable for most oxygen prescriptions.




Cabinet Filter



LCD Display Screen



Flow Meter



Humidifier Bottle

**48 Hour service guarantee with backup machines | 16 service centers across India**

### LOCATIONS

JAMMU | LUDHIANA | CHANDIGARH | DELHI | DEHRADUN | BAREILLY | JAIPUR | LUCKNOW | VARANASI  
AHMEDABAD | KOLKATA | MUMBAI | HYDERABAD | BANGALORE | CHENNAI | KOCHI

Marketed by: Sanrai Med India Pvt. Ltd.  
Office Tower B, Elante Office No.506, Plot No.178-178A,  
Industrial Area Phase 1, Chandigarh, 160002

Website: [www.shop.sanraimed.com](http://www.shop.sanraimed.com)  
Customer Support: +91 7341103371  
Email: [customersupport@sanraimed.com](mailto:customersupport@sanraimed.com)