

Silver Sport 2

TECHNICAL SPECIFICATIONS

Weight Capacity	250 lbs 300 lbs 350 lbs
Height	36"
Seat To Floor Height	17.5" - 19.5"
Carton Shipping Weight	48 lbs 49 lbs 50 lbs
Back Height	16"



Features

- Powder-coated silver vein steel frame (Figure A)
- Embossed vinyl upholstery is durable, lightweight, attractive and easy-to-clean
- Urethane tires mounted on composite wheels provide durability and low maintenance
- 8" front casters, adjustable in three positions
- Padded armrests provide additional comfort
- Carry pocket on backrest provides additional convenience
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Plastic foot plates with heel loops on footrests
- Comes with push-to-lock wheel locks

WE ARE LOCATED

JAMMU LUDHIANA CHANDIGARH DELHI DEHRADUN BAREILLY JAIPUR LUCKNOW VARANASI
AHMEDABAD KOLKATA MUMBAI HYDERABAD BANGALORE CHENNAI KOCHI